

# *In vivo* $\mu$ CT of the myocardium



## *In vivo* quantification of myocardial infarction using ExiTron™ MyoC 8000

Stefan Sawall<sup>1</sup>, Danielle Franke<sup>2</sup>, Anne Kirchherr<sup>2</sup>, Jan Beckendorf<sup>3,4</sup>, Jan Kuntz<sup>1</sup>, Joscha Maier<sup>1</sup>, Alexander Kraupner<sup>2</sup>, Johannes Backs<sup>3,4</sup>, Andreas Briel<sup>2</sup>, and Marc Kachelrieß<sup>1</sup>

<sup>1</sup> Medical Physics in Radiology, German Cancer Research Center (DKFZ), Heidelberg, Germany

<sup>2</sup> nanoPET Pharma GmbH, Berlin, Germany

<sup>3</sup> Department of Molecular Cardiology and Epigenetics, University of Heidelberg, Heidelberg, Germany

<sup>4</sup> German Centre for Cardiovascular Research (DZHK), Heidelberg, Germany

\* Corresponding author (stefan.sawall@dkfz.de)

## Introduction

Cardiovascular disease remains the leading cause of deaths worldwide, accounting for more than 17 million deaths per year<sup>1</sup>. With the exceedingly high number of people affected by this disease, a strong focus is being placed on the advancement of diagnostic imaging techniques that allow for improved detection of cardiovascular pathologies as well as monitoring of potential strategies for disease prevention and treatment. In recent years, *in vivo* imaging in

small animal models of cardiovascular disease plays a central role in cardiovascular research and the models have certainly proven their value and clinical translatability. In this study we use a mouse model of myocardial infarction and develop a  $\mu$ CT method involving the use of the innovative contrast agent ExiTron™ MyoC 8000 to non-invasively detect and effectively quantify myocardial infarction *in vivo*.


## Materials and methods

Experiments were performed on healthy C57BL/6J mice (n = 10) and C57BL/6J mice with myocardial infarction (n = 5) inflicted by surgical ligation of the left anterior descending (LAD) artery<sup>2</sup>. Following surgery, the animals were allowed to recover for a period of 2 weeks before start of the imaging experiments.  $\mu$ CT was performed on all animals after intravenous injection of ExiTron™ MyoC 8000 (Viscover™, nanoPET Pharma GmbH, Berlin, Germany), a nanoparticulate imaging agent having a high iodine content of 210 mg/mL. To investigate contrast agent kinetics, healthy anaesthetized mice were imaged before and after injection of ExiTron MyoC 8000 at a dose of 1050 mg iodine/kg body weight (i.e., 125  $\mu$ L per 25 g mouse) followed by a saline flush (25  $\mu$ L). Imaging was performed repetitively over 250 min with an additional scan after 24 h using a flat detector-based  $\mu$ CT system (Siemens Healthineers, Forchheim, Germany). Each scan was performed over 80 s with a rotation time of 19 s per revolution and image reconstructions were performed using a low-dose phase-correlated reconstruction method<sup>3,4</sup>. ECG data was obtained via electrodes attached to the animal's paw and the respiratory signal was derived using a pneumatic pillow (Small Animal Instruments, Stony Brook, USA). To allow for repetitive image acquisitions, the scan protocol was optimized to minimize radiation dose: each measurement comprised 50 mGy, resulting in an overall radiation dose of 600 mGy. The X-ray tube was operated using a tube voltage of 80 kV and a tube current of 50 mA. Signal enhancement was measured in different regions of interest (ROIs), including the left ventricle (corresponding to vasculature), myocardium, liver, spleen, kidney, brain and muscle.

For detection and subsequent quantification of myocardial infarction, the signal enhancement in the myocardium of infarcted mice was measured at 210 min post injection of the contrast agent. After completion of the imaging experiments, the mice were sacrificed and histological analysis was performed on paraffin-embedded sections of the heart using Masson's Trichrome<sup>5</sup>. Quantification of infarct size in histological stains and reconstructed 3D images was performed as previously described<sup>6</sup>.

## Results and discussion


In order to investigate the kinetics of ExiTron MyoC 8000, signal enhancements in various tissues of healthy mice were measured over a period of 24 h. The time-enhancement curve of the left ventricle (Fig. 1) illustrates that the contrast agent has a long blood half-life of ~2 h in mouse, with the maximum enhancement occurring immediately after injection. At 250 min post injection the signal enhancement in the blood has nearly reached a baseline value denoting that at this time point the contrast agent has practically cleared from the blood. Similar to the signal enhancement in the blood, the signal enhancement in the myocardium rapidly increases following injection of the contrast agent (Fig. 1). The peak enhancement is observed at 20 min post injection but the CT value at this time point is likely overestimated due to an overexposure artifact resulting from the strong contrast in the adjacent ventricles. Thereafter, the signal enhancement in the myocardium remains constant, at least until the 250 min imaging time point. As a result, ExiTron MyoC 8000 enables



**Figure 1:** Time-dependent signal enhancement in the left ventricle and myocardium of healthy mice after injection of ExiTron MyoC 8000 at a dose of 1050 mg iodine/kg body weight.


not only functional cardiac imaging during the first 2 h of the blood pool phase but also imaging of the myocardium at a later phase.  $\mu$ CT images obtained during these two phases highlight the differences in signal enhancement within the various tissues (Fig. 2). The vessels, ventricles and aortic arch can be easily discriminated from surrounding tissues at 15 min post injection, whereas vessel contrast is diminished at 210 min post injection enabling delineation of the myocardium.

Besides providing signal enhancement in the blood and myocardium, ExiTron MyoC 8000 also results in the enhancement of the liver, spleen and kidneys, indicating



**Figure 2:** Axial, coronal and sagittal  $\mu$ CT reconstruction images obtained after injection of ExiTron MyoC 8000 in the healthy mouse at a dose of 1050 mg iodine per kg body weight.


**A.** In the early phase (15 min p.i.), the blood is highly contrasted allowing for functional cardiac imaging (C/W = 500/1500 HU). **B.** In the later phase (210 min p.i.), vessel contrast is diminished enabling delineation of the myocardium (C/W = 500/850 HU).



**Figure 3:** Qualitative comparison of  $\mu$ CT and histology images of mouse with myocardial infarction after injection of ExiTron MyoC 8000. Axial  $\mu$ CT image at 210 min p.i. and corresponding histology image showing myocardial infarction (red arrows). In the histology image, viable tissue is depicted in red and infarct scar as well as local collagenous tissue, in green. Note the viable papillary muscle (red asterisks), which can be identified in both the  $\mu$ CT and histology image.

that these organs play a role in contrast agent elimination<sup>7</sup>. As expected, the muscle and brain tissues do not show any significant uptake of the contrast agent. At the 24 h time point (1440 min) the signal enhancement curves of all measured tissues converge with the baseline indicating that the contrast agent is practically eliminated from the body, offering the possibility of multiple injections<sup>7</sup>.

Since imaging at 210 min post injection of ExiTron MyoC 8000 allows for effective visualization of the myocardium, this time point was selected for subsequent studies in a mouse model of myocardial infarction. Through  $\mu$ CT reconstruction images, infarct sizes were quantified and compared to those obtained via histopathology. A qualitative comparison of the  $\mu$ CT reconstruction images and microscopy images of histological sections depicts that the contrast agent is homogeneously distributed in the healthy myocardium, but is not taken up by areas of myocardial infarction (Fig. 3). The infarct sizes obtained via imaging as well as via histology vary in the range of 14 % and 60 % (Fig. 4), and are similar to those previously reported<sup>6</sup>. A quantitative comparison of infarct sizes obtained via the two methods shows a significantly strong correlation ( $R^2 = 0.98$ ), demonstrating that the developed method allows for the detection and effective quantification of myocardial infarction *in vivo*.



**Figure 4:** Comparison of infarct size obtained via  $\mu$ CT reconstruction images and via histomorphometry shows a significantly strong correlation ( $R^2 = 0.98$ ).

## Conclusion

In this study we develop a  $\mu$ CT method involving ExiTron MyoC 8000 and a low-dose phase-correlated reconstruction algorithm to monitor cardiac processes in small animals *in vivo*. The method allows for functional cardiac imaging as well as imaging of the myocardium, offering an accurate and reproducible method for *in vivo* detection and quantification of infarct size. Thus, the approach presents researchers with an innovative tool for the monitoring of cardiovascular pathologies in small animal models in order to evaluate the response to therapy.

Viscover™ Product	Order No.
ExiTron™ U, 1 x 5 injections	130-095-142
ExiTron™ U, 5 x 5 injections	130-095-143
ExiTron™ V, 1 x 5 injections	130-095-283
ExiTron™ V, 5 x 5 injections	130-095-284
ExiTron™ P, 1 x 5 injections	130-095-144
ExiTron™ P, 5 x 5 injections	130-095-145
ExiTron™ nano 6000, 1 x 5 injections	130-095-146
ExiTron™ nano 6000, 5 x 5 injections	130-095-147
ExiTron™ nano 12000, 1 x 5 injections	130-095-698
ExiTron™ nano 12000, 5 x 5 injections	130-095-700
ExiTron™ MyoC 8000, 1 x 5 injections	130-095-701
ExiTron™ MyoC 8000, 5 x 5 injections	130-095-702

## References

1. World Health Organization (WHO) Fact Sheet, Cardiovascular Diseases (CVDs), May 2017: <http://www.who.int/mediacentre/factsheets/fs317/en/>.
2. Wang, J. et al. (2006) A simple and fast experimental model of myocardial infarction in the mouse. *Tex Heart Inst J.* 33(3): 290-293.
3. Sawall, S. et al. (2011) Low-dose cardio-respiratory phase-correlated cone-beam micro-CT of small animals. *Medical Physics.* 38(3): 1416–1424.
4. Maier, J. et al. (2014) Assessment of dedicated low-dose cardiac micro-CT reconstruction algorithms using the left ventricular volume of small rodents as a performance measure. *Medical Physics.* 41(5): Article ID051908.
5. Zhao, X. et al. (2015) Optical projection tomography permits efficient assessment of infarct volume in the murine heart post myocardial infarction. *American Journal of Physiology-Heart and Circulatory Physiology.* 309(4): H702–H71.
6. Stegger, L. et al. (2006) Accurate noninvasive measurement of infarct size in mice with high resolution PET. *Journal of Nuclear Medicine.* 47(11): 1837–1844.
7. Sawall, S. et al. (2017) In Vivo Quantification of Myocardial Infarction in Mice using Micro-CT and a Novel Blood Pool Agent. *Contrast Media Mol Imaging.* <https://doi.org/10.1155/2017/2617047>.

# Scintica:

[www.scintica.com](http://www.scintica.com)

+1 (519) 914 5495  
[info@scintica.com](mailto:info@scintica.com)

562 Waterloo St. London, ON,  
Canada N6B 2P9