

WEBINAR:

(July 24, 2024) Lab Incubator Decontamination: Why, When and How

Questions and answers from the April 24, 2024 webinar titled: “*Lab Incubator Decontamination: Why, When and How*”

This document includes questions we received and answered during the webinar, as well as those that we did not have time to address.

1. Is there any residual antimicrobial activity after the cycle is complete?

Randall: No, there's no residual antimicrobial activity. Hydrogen peroxide degrades completely and spontaneously too.

2. Does it work for microplasma, and is high hydrogen peroxide? Can it be removed? Can it remove microplasma?

Randall: Yes, it is effective against most mycoplasma species.

3. How often should users decontaminate their incubator?

Randall: We recommend once a month, depending on how many incubators you have in your lab, and what the workload is. But I would say, start with once a month.